



Thalassa S.E. Valkenburg MSc, MSc (tsev1@le.ac.uk)¹, Dr R.L. Cordell¹, Prof P.S. Monks¹ and Dr H.C. Pandya²

1. Department of Chemistry, University of Leicester, 2. Department of Immunity, Infection & Inflammation, University of Leicester

Salbutamol Screening

Salbutamol is a medication used to relieve the symptoms of asthma, bronchospasm and chronic obstructive pulmonary disease (COPD).




- First line therapy: inhalation.
- Next step: intravenous administration.

At present there are insufficient data to recommend a dosage regime for children and extrapolation of data from one age group to another is recognized to be unsafe¹. Excessive administration can lead to lactic acidosis, which precipitates respiratory failure².

- Research gap: dose finding studies to determine optimal doses for paediatric asthma applications.
- Challenge: measuring blood concentrations of a drug in children.
- Research aim: development of a novel assay to detect and quantify salbutamol from micro-volumes of blood.

Blood Samples

- Dried Blood Spots (DBS): spotting and drying of micro-volumes blood onto filter paper³.
- Sampling issues using filter paper: haematocrit variation, sample inhomogeneity and unknown sample volumes.



Other blood collection devices have been introduced. One of them is used alongside filter paper in this study; Volumetric Absorptive MicroSampling (VAMS)*, also named Mitra⁴.

*kindly donated by James B. Rudge, Phenomonex

- Research samples: 10 µl blood samples spiked with salbutamol directly after spotting and with d3-salbutamol as internal standard after drying and extraction. The samples were evaporated to dryness and derivatised with BSTFA plus 1% TMCS prior to analysis with a developed GC-MS method.

Assay Performance

- Lower limit of quantification (LLOQ): 5 ng/ml - below the therapeutic range.
- Accuracy and precision: <20% for LOD samples and <15% for other samples, meeting the guidelines for bioanalytical method validation of the United States Food and Drug Administration (US FDA)⁵.
- Recovery: salbutamol recovered from DBS shows a higher recovery (86%) compared to blood sampled with the Mitra device (68%). This applies, however, only to samples left to dry for 2 hours and processed immediately upon drying.
- Stability: DBS samples show a decreasing recovery with increasing drying time compared to Mitra tip samples (Fig 1).
- Application test: three healthy volunteers were administered 1 mg of salbutamol via a 100 µg per dose inhaler and spacer (Fig 2).




Figure 1: Effect on salbutamol recovery from DBS samples of increasing drying/storage times.




Figure 2: Screening of salbutamol over time from 10 µl blood sampled with Mitra tips.

Conclusion

- An assay to detect and quantify salbutamol from micro-volumes of blood has been developed, optimised and validated according to US FDA guidelines, especially valuable for paediatric asthma applications.
- Two blood collection devices were tested and compared for this study; the novel developed VAMS devices have demonstrated sampling and stability advantages compared to the conventional filter cards used for DBS research and applications.
- As a final method evaluation, salbutamol could be detected in samples from three healthy volunteers dosed with salbutamol.

References

1. Cella M., Knibbe C., Danhof M. and Della Pasqua O. (2010). What is the right dose for children? *British Journal of Clinical Pharmacology*.
2. Tobin A. (2005). Intravenous salbutamol: too much of a good thing? *Critical Care and Resuscitation*.
3. Mei J.V., Alexander R., Adam B.W. and Hannon W.H. (2001). Use of filter paper for the collection and analysis of human whole blood specimens. *The Journal of Nutrition*.
4. Denniff P. and Spooner N. (2014). Volumetric absorptive microsampling: a dried sample collection technique for quantitative bioanalysis. *Analytical Chemistry*.
5. US Department of Health and Human Services, Food and Drug Administration, Center for Drug Evaluation and Research, Center for Veterinary Medicine (2011). Guidance for Industry, Bioanalytical Method Validation, <http://www.fda.gov/downloads/Drugs/Guidances/ucm070107.pdf>.

